

MISHIN GENERATOR WITH CAPACITIVE COILS

MANUAL

The complete set (capacitive coils) includes:



The complete set of Mishin Coil (Mishin generator)

- Sine autogenerator.
- High quality power supply 5V 2000mA, with output filtering, cable (length is 100cm).
- Toroidal capacitive coil with a diameter of 6.5cm. TMA (Torus of Mishin Alexandr) is a round coil shaped like a tor. It is smaller than the radius of the zone of influence, but greater than the density of the impact, in contrast to the flat coil, which has a smaller impact density and greater radius of influence. It has a stronger field than DMA, so you should not use it for a long time during one session. The best option is to use DMA then TMA, then again DMA and then TMA and so on. The coil operates at 300 kHz
- Disk capacitive coil of wires in varnish 0.5mm. DMA (Mishin Alexander Disk) is a wooden round disc (flat coil). The disc on both sides are covered with plywood 4mm,

which when touched to the body does not change and does not distort the frequency of the coil. The length of the cable from the coil to the generator is 2 m. The coil operates at 300 kHz.

– Static field indicator.

– USB-cable that allows the generator to work from any power unit with USB (a power bank, a laptop, etc.)

Getting the Device Ready for the Treatment Session

- 1) Connect the coil (a wooden round disk) into the sine autogenerator. You should plug the coil into a red socket (on the picture it is named “ВЫХОД СИГНАЛА”).



Picture 1. the sine autogenerator


- 2) Connect the device power supply with the sine autogenerator using the socket for power connection shown on the Picture 2



Picture 2

- 3) Plug the device into an outlet or into a power bank or similar device using USB-cable. When you do that two LED indicators will appear on the panel of the

device: the red LED indicates the presence of power, the green LED lights up when the coil is connected and indicates that the device is ready for being used.



Picture 3

Treatment Sessions:

- 1) Remove all metal products (chains, rings, etc.) and put all metal products away. (especially closed jewelry such as chains, bracelets, etc. If a necklace is worn on the neck, the coil action will not reach the head.
- 2) Before the session, drink one glass of water.
- 3) Put a coil on the organ you are treating. You can apply the coils by any side. To increase the effectiveness of a torus or disk on the top of the coil, you can put your hand. All diseases have a closed ring structure. The difference is only in size. Very coarse coils will work with very small viruses, so we use a torus. There should be at least two coils. This is a flat disk (medium or large) and a torus. With the full coincidence of the phases (voltage and current) of the signal supplied by the generator to the coil, we physically get a standing wave, when the phases shift, it starts to rotate like a "fan". Up to 90 degrees of phase shift implosion processes occur in the coils, which benefit us. It also occurs with an increase in the frequency, higher than the operating frequency of the coil (the processes in the ES harmonics are implosion), and when the frequency of the harmonics decreases, they emit radiation, and the radiation around us is already full. At 90 degrees, the rotation speed will be maximum and the flow merges into a ring.
- 4) Time of the session depends on the disease you are treating. Overdoing in treatment is not necessary. Sessions should last up to 15-60 minutes a day or every other day (according to the sensations). To help the kidneys remove toxins (and they are usually clogged in all) we put on each kidney a torus every day for 10 minutes. Up to 60-90 minutes of the coil (torus or disk) operation per day is

sufficient for the sore removal process. The main thing is not to "overdose", so that the body has time to regenerate and remove toxins.

The speed of recovery depends on several factors:

1. Resources of the body: how much illnesses have worn out the life support system.
2. Locations of a patient. We live in an environment with constant pulsation from electrical devices: mobile phones, tablets, computers, cell towers. These pulsations are the cause of many diseases. Therefore, the process of recovery is better to start outside the city, where the sources of such pulsations are less.

Note! Here is a table of the diseases and time, it is not complete, so if your disease is not listed here, please, contact us and we will provide you with more information. The time given in the table is not advisable but not strict, which means you can use it a bit longer like not 15 but 20 min or so.

List of Diseases and Recommended Treatment Sessions

COLD AND FLU

Flat coil or torus – 20 minutes on the chest and 20 minutes on the kidneys. Do it for 2-3 days in the morning and in the evening. Advice: Drink acidified teas.

ANGINA

Torus – 15-20 minutes at the throat. Do it for 2-3 days in the morning and in the evening.

PHYSICAL INJURY

Flat coil (can be alternated with a torus) – 15-20 min on the injured place. Use it 1-2 times a day to healing.

ARTROZ, OSTEOARTROZ, OSTEOPOROZ, ARTRITY, ARTHROSIS

These are cholesterol plaques that grow, often, from insufficient blood supply to the

joints.

It all depends on the specific case. Flat coil – 30-40 minutes or torus – 10 min on the chest and stomach.

Do it for 2-3 months a day every other day. After 2 weeks of use, take one week break.

PROSTATITIS, PROSTATE ADENOMA

The cause of prostatitis and similar problems is the problem of filtration of blood by the kidneys. Because of this, toxins accumulate.

Accumulation of toxins causes irritation of the prostate. It, like a sponge, absorbs it all.

Flat coil (can alternate with the torus) – 30-40 min on the kidneys, the lower abdomen and coccyx. Do it in a day to the complete recovery.

DIABETES MELLITUS I AND II

Flat coil – 30-40 minutes torus – 10 min on the chest and stomach. Do it in a day. After 2-3 weeks of use, take a week break.

KIDNEYS, KIDNEYSTONES, PYELONEPHRITIS

Back pain comes from the kidneys. Flat coil – 20-30 minutes on the lower back or torus – 10 min on each kidney.

Do it 1-2 times a day it in a day. The color and smell of urine indicates the number of toxins in the blood. Take breaks for a week with enough darkening.

GALLSTONES, GALLSTONES, CHOLECYSTITIS

Torus – 15 min in the area of the gallbladder. Do it in a day. After 2-3 weeks of use, take one week break.

ALLERGIES

Flat coil – 30-40 minutes on the chest and stomach. Do it in a day. After 2 weeks of use, take one week break.

VARICOSE VEINS

Flat coil – 30 min on problematic places. Use it several times a day to complete recovery.

SCARS, SCARTISSUE

Flat coil (can be alternated with the torus) – 15-20 minutes on the scared areas. Use it 1-2 times a day to complete recovery.

HEADACHES, MIGRAINES

Many problems with headaches are associated with improper blood supply. If the headaches do not go away: Torus – 10 minutes on every kidney. Use it once a day.

NICOTINE ADDICTION

Flat coil – 30-40 minutes on the back. Use it once a day. Toxic removes in an hour, can pierce a cough.

POLYCYSTIC OVARY SYNDROME (PCOS)

Your treatment plan with Mishin coil should be: on the first day use torus on your kidneys for 5 min each, then put it on your liver for 10 min; the next day use a flat coil for 10-15 min on places with ovarian cysts. Do it for a week. After the first week you should put a flat coil for 30-40 minutes on problematic places in a day. Mishin coil is a perfect natural remedies for polycystic ovary syndrome

HERNIAS, PROTRUSIONS OF THE VERTEBRAL DISCS

Pain must be removed, but we cannot restore mechanics when it is directly deformed. Try a flat coil (you can alternate with a torus) 30-40 minutes in the problem areas. Use it in a day.

VSD, PANIC ATTACKS

Flat coil – 30-40 minutes torus – 10 min on the chest, abdomen. Use it in a day.

DISEASES OF THE SKIN – VITILIGO, WHEN WHITE SPOTS APPEAR ON THE SKIN

Flat coil – 30-40 minutes on the chest, abdomen. Try additionally for 10-15 minutes on problem areas.

Use it once a day.

JAUNDICE

Flat coil – 30-40 minutes on the chest, stomach, liver. Use it in a day.

PAPILLOMAS, POLYPS, WARTS (COLORLESS BIRTHMARKS AND SKIN TAGS)

Flat coil – 30-40 min on the chest, abdomen or torus – 15 min on problematic areas.

Use it in a day to complete recovery

ULCERS OF THE STOMACH AND DUODENUM

In the morning, at lunch and in the evening for half an hour before meals take a glass of water (not cold, room temperature).

Dissolve in it three teaspoons of dry potato starch and add five six droplets of apothecary 5% iodine. Stir everything and drink.

Use a glass three times a day in a week. Flat coil – 30-40 minutes on the chest, abdomen. Use it in a day

ANGINA PECTORIS AND ISCHEMIC HEART DISEASE. COMPLICATION OF THE HEART AFTER THE FLU

Flat coil – 30-40 minutes or torus – 10 min on the chest, abdomen. Use it in a day

MYOCARDIAL INFARCTION

Flat coil for one week – 30-40 minutes on the chest, stomach in a day.

Then Tor for one week – 30-40 minutes on the chest, stomach in a day.

At first there will be a pulling pain.

DISEASES OF THE THYROID GLAND, NODES ON THE THYROID GLAND

Torus (can be alternated with flat) – 15-20 min in the thyroid area. Use it in a day.

HEPATITIS AND AIDS

Try a flat coil – 30-40 minutes on the chest, stomach in a day. There are no statistic data on the results of AIDS treatment yet.

OVARIANCYST

Flat coil (can be alternated with the torus) – 30-40 min on problematic places. Use it in a day

FUNGAL DISEASES

Flat coil (can be alternated with the torus) – 20 min on problematic places. Use it in a day

TEETH, GUMS, PERIODONTAL DISEASE

Flat coil (can be alternated with the torus) – 15-20 minutes in the mouth area. Use it 1-2 times a day to complete recovery.

HEART ARYTHMY

Flat coil – 30-40 minutes on the chest, abdomen. Use it in a day

CANCER

The first stages: Flat coil- 30-40 minutes torus – 10 min on the chest, abdomen. Use it in a day.

There are stages when the defeat of the internal organs has passed a critical point.

In these cases, the device can only reduce suffering.

HERNIATED HERNIA

Flat coil – 30-40 minutes on the chest. Use it for 2-3 months in a day.

DRUGD ADDICTION

Flat coil – 30-40 minutes on the chest, abdomen.

If someone has used drugs, the device removes the dependence of toxic substances.

PSORIASIS

Flat coil – 30-40 minutes on the chest, abdomen. Use it in a day.

LYMPHOSTASIS

Flat coil- 30-40 minutes on the chest, abdomen. Use it for 5-10 days in a day

ASTHMA OF ALLERGIC ORIGIN

Flat coil – 30-40 minutes on the chest and on the back. Use it for 3 months in a day.

ISCHEMIC STROKE OF THE BRAIN

Flat coil – on the chest and stomach for 30-40 min + flat coil on the head area for 5-10 min. Use it in a day.

CYSTOPANCREAS

Flat coil – 30-40 minutes on the chest, abdomen. Use it for 2-3 weeks in a day.

CATARACT, MACULAR DEGENERATION OF THE RETINA

Most often there is an improvement within the first week, but it takes time to recover completely.

Try a flat coil – 5-10 minutes on the head.

Use it once a day. Dystrophy of the retina: the recovery processes are slowly going on and for the time being, there are only a few cases where a person began to see better after a brief impact.

PARKINSON

Flat coil – 30-40 minutes on the chest, stomach + 5-10 minutes on the head. Use it in a day.

MULTIPLESCLEROSIS

Flat coil – 30-40 minutes on the chest, stomach. Use it in a day.

EPILEPSY

Flat coil – 30-40 minutes on the chest, stomach. Use it in a day.

DISEASES OF THE BLOOD. IN PARTICULAR, LYMPHATIC LEUKEMIA

Flat coil – 30-40 minutes on the chest and on the back. Use 5-7 consecutive days, 1 time per day.

Pay attention to the sensations in the kidneys and the color of urine.

TOURETTE'S SYNDROME

Flat coil – 30-40 minutes on the chest, stomach. Use it in a day

MULTIPLE SCLEROSIS AND SCLERODERMA

Flat coil – 30-40 minutes on the chest, stomach. Use it in a day

THROMBOSIS OF THE VENOUS SINUS OF THE BRAIN

Flat coil – 30-40 minutes on the chest, stomach + 5-10 minutes on the head. Use it in a day

NEURALGIA OF THE FACIAL TRIGEMINAL NERVE

There was one patient with such a problem. She had great pains, she could not even sleep.

In addition, she had metastases and 2 chemotherapies behind. She was all shaking and could not even speak normally. After 3 sessions for an hour with a flat coil, the main

neurology began to decrease, and pain from the ternary decreased significantly after 20-30 minutes of a flat coil.