

MISHIN GENERATOR WITH A STATIC COIL

MANUAL

The complete set (static coil or Tesla) includes:



- Sine autogenerator.
- High quality power supply 5V 2000mA, with output filtering, cable (length is 100cm).
- A cable for connection of an autogenerator to Power Bank.
- The disk static coil from a wire in a varnish of 0.4mm.

– Static field indicator.

– USB-cable that allows the generator to work from any power unit with USB (a power bank, a laptop, etc.)

The discs on both sides are covered with plywood 4mm, which when touched to the body does not change and does not distort the frequency of the coils.

The length of the cable from the coils to the generator is 2 m.

For new coils, only a sinusoidal signal is used at the output of the autogenerator.

Getting the Device Ready for the Treatment Session

- 1) Connect the coil (a wooden round disk) into the sine autogenerator. You should plug the coil into a red socket (on the picture it is named “ВЫХОД СИГНАЛА”).



Picture 1. the sine autogenerator


- 2) Connect the device power supply with the sine autogenerator using the socket for power connection shown on the Picture 2



Picture 2

- 3) Plug the device into an outlet or into a power bank or similar device using USB-cable. When you do that two LED indicators will appear on the panel of the

device: the red LED indicates the presence of power, the green LED lights up when the coil is connected and indicates that the device is ready for being used.



Picture 3

Treatment Sessions:

- 1) Remove all metal products (chains, rings, etc.) and put all metal products away. (especially closed jewelry such as chains, bracelets, etc. If a necklace is worn on the neck, the coil action will not reach the head.
- 2) Before the session, drink one glass of water.
- 3) Put a coil on the organ you are treating. You can apply the coil by any side. Time of the session depends on the disease you are treating. Overdoing in treatment is not necessary. Sessions should last up to 15-60 minutes a day or every other day (according to the sensations). Up to 60-90 minutes of the coil operation per day is sufficient for the sore removal process. The main thing is not to "overdose", so that the body has time to regenerate and remove toxins.

The speed of recovery depends on several factors:

1. Resources of the body: how much illnesses have worn out the life support system.
2. Locations of a patient. We live in an environment with constant pulsation from electrical devices: mobile phones, tablets, computers, cell towers. These pulsations are the cause of many diseases. Therefore, the process of recovery is better to start outside the city, where the sources of such pulsations are less.

Note! Here is a table of the diseases and time, it is not complete, so if your disease is not listed here, please, contact us and we will provide you with more information. The

time given in the table is not advisable but not strict, which means you can use it a bit longer like not 15 but 20 min or so.

List of Diseases and Recommended Treatment Sessions

COLD AND FLU

Static coil – 20 minutes on the chest and 20 minutes on the kidneys. Do it for 2-3 days in the morning and in the evening. Advice: Drink acidified teas.

ANGINA

Static coil – 20-30 minutes at the throat. Do it for 2-3 days in the morning and in the evening.

PHYSICAL INJURY

Static coil – 15-20 min on the injured place. Use it 1-2 times a day to healing.

ARTROZ, OSTEOARTROZ, OSTEOPOROZ, ARTRITY, ARTHROSIS

These are cholesterol plaques that grow, often, from insufficient blood supply to the joints.

It all depends on the specific case. Static coil – 30-40 minutes on the chest and stomach.

Do it for 2-3 months a day every other day. After 2 weeks of use, take one week break.

PROSTATITIS, PROSTATE ADENOMA

The cause of prostatitis and similar problems is the problem of filtration of blood by the kidneys. Because of this, toxins accumulate.

Accumulation of toxins causes irritation of the prostate. It, like a sponge, absorbs it all.

Static coil – 30-40 min on the kidneys, the lower abdomen and coccyx. Do it in a day to the complete recovery.

DIABETES MELLITUS I AND II

Static coil – 30-40 minutes on the chest and stomach. Do it in a day. After 2-3 weeks of use, take a week break.

KIDNEYS, KIDNEYSTONES, PYELONEPHRITIS

Back pain comes from the kidneys. Static coil – 20-30 minutes on the lower back.

Do it 1-2 times a day it in a day. The color and smell of urine indicates the number of toxins in the blood. Take breaks for a week with enough darkening.

GALLSTONES, GALLSTONES, CHOLECYSTITIS

Static coil – 30 min in the area of the gallbladder. Do it in a day. After 2-3 weeks of use, take one week break.

ALLERGIES

Static coil – 30-40 minutes on the chest and stomach. Do it in a day. After 2 weeks of use, take one week break.

VARICOSE VEINS

Static coil – 30 min on problematic places. Use it several times a day to complete recovery.

SCARS, SCARTISSUE

Static coil – 15-20 minutes on the scared areas. Use it 1-2 times a day to complete recovery.

HEADACHES, MIGRAINES

Many problems with headaches are associated with improper blood supply. Use it once a day.

NICOTINE ADDICTION

Static coil – 30-40 minutes on the back. Use it once a day. Toxic removes in an hour, can pierce a cough.

POLYCYSTIC OVARY SYNDROME (PCOS)

Your treatment plan with Mishin coil should be: on the first day use a static coil on your kidneys for 20 min each, then put it on your liver for 10 min; the next day use a static coil for 10-15 min on places with ovarian cysts. Do it for a week. After the first week you should put a static coil for 30-40 minutes on problematic places in a day. Mishin coil is a perfect natural remedy for polycystic ovary syndrome

HERNIAS, PROTRUSIONS OF THE VERTEBRAL DISCS

Pain must be removed, but we cannot restore mechanics when it is directly deformed. Try a static coil 30-40 minutes in the problem areas. Use it in a day.

VSD, PANIC ATTACKS

Static coil – 30-40 minutes on the chest, abdomen. Use it in a day.

DISEASES OF THE SKIN – VITILIGO, WHEN WHITE SPOTS APPEAR ON THE SKIN

Static coil – 30-40 minutes on the chest, abdomen. Try additionally for 10-15 minutes on problem areas.

Use it once a day.

JAUNDICE

Static coil – 30-40 minutes on the chest, stomach, liver. Use it in a day.

PAPILLOMAS, POLYPS, WARTS (COLORLESS BIRTHMARKS AND SKIN TAGS)

Static coil – 30-40 min on the chest, abdomen. Use it in a day to complete recovery

ULCERS OF THE STOMACH AND DUODENUM

In the morning, at lunch and in the evening for half an hour before meals take a glass of water (not cold, room temperature).

Dissolve in it three teaspoons of dry potato starch and add five six droplets of apothecary 5% iodine. Stir everything and drink.

Use a glass three times a day in a week. Static coil – 30-40 minutes on the chest, abdomen. Use it in a day

ANGINA PECTORIS AND ISCHEMIC HEART DISEASE. COMPLICATION OF THE HEART AFTER THE FLU

Static coil – 30-40 minutes on the chest, abdomen. Use it in a day

MYOCARDIAL INFARCTION

Static coil for one week – 30-40 minutes on the chest, stomach in a day. At first there will be a pulling pain.

DISEASES OF THE THYROID GLAND, NODES ON THE THYROID GLAND

Static coil – 30-40 min in the thyroid area. Use it in a day.

HEPATITIS AND AIDS

Try a static coil – 30-40 minutes on the chest, stomach in a day. There are no statistic data on the results of AIDS treatment yet.

OVARIANCYST

Static coil – 30-40 min on problematic places. Use it in a day

FUNGAL DISEASES

Static coil – 20 min on problematic places. Use it in a day

TEETH, GUMS, PERIODONTAL DISEASE

Static coil – 15-20 minutes in the mouth area. Use it 1-2 times a day to complete recovery.

HEART ARYTHMY

Static coil – 30-40 minutes on the chest, abdomen. Use it in a day

CANCER

The first stages: Static coil - 30-40 minutes on the chest, abdomen. Use it in a day.

There are stages when the defeat of the internal organs has passed a critical point.

In these cases, the device can only reduce suffering.

HERNIATED HERNIA

Static coil – 30-40 minutes on the chest. Use it for 2-3 months in a day.

DRUGD ADDICTION

Static coil – 30-40 minutes on the chest, abdomen.

If someone has used drugs, the device removes the dependence of toxic substances.

PSORIASIS

Static coil – 30-40 minutes on the chest, abdomen. Use it in a day.

LYMPHOSTASIS

Static coil - 30-40 minutes on the chest, abdomen. Use it for 5-10 days in a day

ASTHMA OF ALLERGIC ORIGIN

Static coil – 30-40 minutes on the chest and on the back. Use it for 3 months in a day.

ISCHEMIC STROKE OF THE BRAIN

Static coil – on the chest and stomach for 30-40 min + static coil on the head area for 5-10 min. Use it in a day.

CYSTOPANCREAS

Static coil – 30-40 minutes on the chest, abdomen. Use it for 2-3 weeks in a day.

CATARACT, MACULAR DEGENERATION OF THE RETINA

Most often there is an improvement within the first week, but it takes time to recover completely.

Try a static coil – 5-10 minutes on the head.

Use it once a day. Dystrophy of the retina: the recovery processes are slowly going on and for the time being, there are only a few cases where a person began to see better after a brief impact.

PARKINSON

Static coil – 30-40 minutes on the chest, stomach + 5-10 minutes on the head. Use it in a day.

MULTIPLESCLEROSIS

Static coil – 30-40 minutes on the chest, stomach. Use it in a day.

EPILEPSY

Static coil – 30-40 minutes on the chest, stomach. Use it in a day.

DISEASES OF THE BLOOD. IN PARTICULAR, LYMPHATIC LEUKEMIA

Static coil – 30-40 minutes on the chest and on the back. Use 5-7 consecutive days, 1 time per day. Pay attention to the sensations in the kidneys and the color of urine.

TOURETTE'S SYNDROME

Static coil – 30-40 minutes on the chest, stomach. Use it in a day

MULTIPLE SCLEROSIS AND SCLERODERMA

Static coil – 30-40 minutes on the chest, stomach. Use it in a day

THROMBOSIS OF THE VENOUS SINUS OF THE BRAIN

Static coil – 30-40 minutes on the chest, stomach + 5-10 minutes on the head. Use it in a day

NEURALGIA OF THE FACIAL TRIGEMINAL NERVE

There was one patient with such a problem. She had great pains, she could not even sleep.

In addition, she had metastases and 2 chemotherapies behind. She was all shaking and could not even speak normally. After 3 sessions for an hour with a static coil, the main neurology began to decrease, and pain from the ternary decreased significantly after 20-30 minutes of a static coil.